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LED Controller PLC300 – Standard

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Operating Instructions

LED Controller PLC300 – Standard

- 1 Safety instructions
- 2 Usage and functions
- 3 Connection options
- 4 Operation
- 4.1 Display & controls
- 4.2 Modes & functions
- 5 Troubleshooting
- 6 Environment and disposal
- 7 Warranty
- 8 Liability

IMPORTANT!

Read carefully before use! Store for later reference!

The content of this document has been checked for compliance with the products described. Nevertheless, deviations cannot be excluded, so that no guarantee can be given for full compliance. The manufacturer is not liable for damages due to missing or incorrect information. Changes or additions are always implemented without notice in the following issues.

1 Safety instructions

The instructions in all manuals must be followed. Failure to follow this instruction may result in property damage and personal injury or even mortal danger. With this controller, only the designated LED lights may be operated. (See → 2)
Switch off the LED controller before changing the pin assignment.

Never look directly into the LED bulbs when switching on or during operation!

2 Usage and functions

The controller is used to connect and control Promicron LED spots, LED flashlights and LED ring lights.

3 Connection Options


The following table informs you about the possible combinations of LED bulbs. If both outputs (A and B) of channel 1 are equipped with LED spots, these are operated synchronously. An LED spot and another LED illuminant cannot be operated simultaneously on one channel.

'X' = connection possible / '-' = assignment not compatible

* = two flashlights with Y adapter conY-2linePlug synchronously controllable

** = Alternatively, flashlights with Y-adapter conY-2lineMiniDIN connectable

Assignment options	channel:	1 A	1 B	2
coolSpot (all versions)	plug 3,5"	X	-	-
coolLine illuminator *	plug 3,5"	X	-	X**
coolRing 66mm, 40 LED	Mini-DIN 9pol	-	-	X

4 Operation

4.1 Display & controls

see photo

4.2 Modes & functions

Pushbutton / rotary knob **channel 1**
press : LED Controller ON / OFF
turn/rotate : brightness channel 1

Pushbutton / rotary knob **channel 2**
turn/rotate : brightness channel 2
multiple assignment pushbutton (LED flashes yellow)
1x press : select segment pattern with rotary knob
2x press : rotate segment pattern with rotary knob

Special functions pushbutton channel 2

(yellow LED is blinking constantly)
Pressing the button once, activates the segment pattern selection. Turning the rotary knob sets six preset segment patterns one after the other. Press the button again to rotate the selected segment pattern with the knob. Press again to switch on the brightness control via the rotary knob. In order to skip one of the three functions, the pushbutton can also be pressed twice in succession. In the default position (brightness setting) the yellow LED does not flash.

5 Troubleshooting

Error : the LED lights are not working
Cause : spot defect or controller defective
Fix : please check spot on other channel; replace if necessary

6 Environment and disposal

The crossed-out wheeled bin symbol means that the product must be disposed of separately from household waste.

7 Warranty

We provide warranty within the scope of the legal provisions. Please return the device in the event of a fault together with a detailed description of the fault and return it to us with sufficient franking.

8 Liability

The product has been designed and manufactured for the intended use described in this document. Any other use is considered improper and may result in damage to the product or personal injury for which the manufacturer cannot be held responsible.

Manufacturer : Promicron AG, Bachmühlweg 24
74366 Kirchheim/Neckar, Germany

Instructions : VS202601EN